

POLY BIOHAZARD WASTE CANS

Safe collection and segregation of potential biohazards

- Ideal for temporary storage of biohazard materials to reduce exposure to hazardous pathogens
- Foot pedal and self-closing lid provide hands-free operation to reduce contamination
- Heavy-wall HDPE construction won't rust or dent and is resistant to most chemical compounds for increased durability
- No seam, leak proof design with molded handles for easy handling
- Flared base provides added stability
- FM approved and meets OSHA and NFPA Code 30 requirements for use with alcohols and other flammable liquids
- Includes OSHA biohazard label on front; use optional label kit 25880 for additional labels to add to the sides, to meet California Health and Safety Code



Compliance:

OSHA 29 CFR 1910.1030



With US and Imported Parts

Capacity	Dimensions W x D x H	Approvals / Compliance	Color	Net Weight	Models
6 gal (22.8 L)	16.5" x 13.5" x 16" (42 x 34.3 x 40.7 cm)	FM, OSHA, NFPA 30	Red	9 lb (4.1 kg)	943BIO
10 gal (37.9 L)	22" x 18" x 18" (55.9 x 45.8 x 45.8 cm)			14 lb (6.4 kg)	945BIO
14 gal (53 L)	22" x 18" x 21" (55.9 x 45.8 x 53.4 cm)			16 lb (7.3 kg)	947BIO

Accessories

Description	Models
Biohazard label kit; for California Health and Safety Code compliance. Includes three labels and instructions.	25880

