



Piggyback Corrosive Cabinet

Starting at
\$995.00

Product Images



Short Description

- Expand storage and safely segregate corrosives
- Feature exclusively formulated paint and poly trays that resist acids, bases and solvents
- Offer same quality safety features as cabinets for flammables
- Available in 17- or 12-gallon capacities
- ChemCor coating options available for harsher corrosive chemicals
- FM-approved
- Ten-year warranty

Description

Research environments using many different types of chemicals have special storage requirements. Whether using acids, bases or solvents, it is critical that incompatible liquids are segregated to avoid disastrous explosions or fire. The SDS (Safety Data Sheet) should be reviewed to determine chemical characteristics and recommended storage practices.

Designed to mount on top of our standard 30- and 45-gallon (114- and 170-L) safety cabinets, these Piggyback cabinets create a safe storage space without taking up any valuable floor space. Cabinets can also be used independently as a wall-mount unit or on a countertop.

Justrite blue steel Piggyback safety cabinets for corrosives have the same quality safety features as the cabinets for flammables. All models come with double-wall construction, dual vents, grounding wire connections, adjustable shelves, leakproof sills, three-point self-latching doors and leveling feet. Self-close door(s) shuts and latches automatically when a fusible link melts at 165°F (74°C) under fire conditions. Unique, concealed self-closing mechanism offers obstruction-free access to contents.

Additionally, they include polyethylene trays attached to galvanized steel shelves and a separate polyethylene liner for the bottom sump resist aggressive chemicals. Easily remove the liner for quick cleaning of drips and leaks. An all epoxy baked-on powder-coat finish, inside and out, provides increased chemical resistance.

For environments using harsher corrosive chemicals, the ChemCor® thermoplastic coating on all interior surfaces provides superior protection against corrosion from drips and off-gassing from harsh chemicals. Seam-free surface eliminates the need for rivets or fasteners which risk corrosion from seep-through chemicals. Non-porous surface resists staining and cleans up easily.

Corrosive Piggyback cabinets meet OSHA and NFPA codes. Cabinets are FM approved and come with a ten-year warranty.

Specifications

Model No	89021-PG
----------	----------

Prop 65

⚠️WARNING:This product can expose you to chemicals including lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Product Options

Door Style:	Self-close
	Manual
Capacity:	12 gallon
	17 gallon
Height:	18
	24
ChemCor:	Yes
	No

