



Slimline Corrosive Cabinet

Starting at
\$1,703.00

Product Images





Short Description

- Ideal for areas with limited floor space or to expand current storage conditions
- Offer same quality safety features as cabinets for flammables
- FM-approved
- Ten-year warranty
- ChemCor coating options available for harsher corrosive chemicals

Description

Research environments using many different types of chemicals have special storage requirements. Whether using acids, bases or solvents, it is critical that incompatible liquids are segregated to avoid disastrous explosions or fire. The SDS (Safety Data Sheet) should be reviewed to determine chemical characteristics and recommended storage practices.

Use Justrite blue Slimline-style cabinets for corrosives to expand storage capacity; its modular size fits nicely alongside existing cabinets. These cabinets have the same quality safety features as the cabinets for

flammables. All models come with double-wall construction, dual vents, grounding wire connections, adjustable shelves, leakproof sills, three-point self-latching doors and leveling feet. Self-close door(s) shuts and latches automatically when a fusible link melts at 165°F (74°C) under fire conditions. Unique, concealed self-closing mechanism offers obstruction-free access to contents.

Additionally, they include polyethylene trays attached to galvanized steel shelves and a separate polyethylene liner for the bottom sump resist aggressive chemicals. Easily remove the liner for quick cleaning of drips and leaks. An all epoxy baked-on powder-coat finish, inside and out, provides increased chemical resistance.

For environments using harsher corrosive chemicals, the ChemCor® thermoplastic coating on all interior surfaces provides superior protection against corrosion from drips and off-gassing from harsh chemicals. Seam-free surface eliminates the need for rivets or fasteners which risk corrosion from seep-through chemicals. Non-porous surface resists staining and cleans up easily.

Slimline cabinets are available in 22- and 54-gallon capacities. They meet OSHA and NFPA codes, are FM approved and come with a ten-year warranty.

Specifications

Model No	89021-SL
----------	----------

Prop 65

⚠️WARNING:This product can expose you to chemicals including lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Product Options

Door Style:	Self-close
	Manual
Capacity:	22 gallon
	54 gallon
ChemCor:	Yes
	No

