



Under Counter Flammable Safety Cabinet

Starting at
\$1,161.00

Product Images





Short Description

- **Efficient Design:** Fits under countertops near point of use, with a recessed toe kick for enhanced ergonomics
- **Durable Construction:** Built with all-welded 18-gauge double-wall steel for superior fire resistance
- **Spill Management:** Adjustable shelves direct spills to a leakproof sump
- **Secure Closure:** Three-point stainless steel bullet latch doors with U-Loc™ handles, available in manual or self-closing options
- **Compliance and Warranty:** Meets NFPA and OSHA regulations, FM-approved, with a ten-year warranty
- **Configurability:** Available in five colors (yellow, red, gray, White, neutral)

Description

Eagle Undercounter Safety Cabinets provide a convenient and safe way to store hazardous liquids near work areas, utilizing undercounter space efficiently. These cabinets feature a recessed toe kick, allowing users to stand closer to the work surface, enhancing ergonomics while mounting flush with other cabinetry. Each cabinet has a capacity of 22 gallons and measures 35 inches in width, 35 inches in depth, and 22 inches in height, making them ideal for fitting under countertops and work benches.

Constructed with all-welded 18-gauge double-wall steel and 1.5 inches of insulating air space, these cabinets offer superior fire resistance and compliance with NFPA and OSHA regulations. The cabinets include adjustable shelves that direct spills to a leakproof 2-inch sump, ensuring safety and compliance. They also feature three-point stainless steel bullet latch doors and U-Loc™ handles for added security. Self-leveling feet adjust to accommodate uneven floors, ensuring stability and alignment.

To configure your undercounter safety cabinet, select from various color options, including yellow for flammables, red for combustibles, gray and White for waste storage, or brown/neutral for lab environments. Choose between manual closing or self-closing doors to meet your specific safety and compliance needs.

Self-closing doors are a valuable feature in safety cabinets, offering enhanced safety and compliance. These doors automatically close and latch in the event of a fire or emergency, providing a critical barrier against the spread of fire and reducing the risk of ignition outside the cabinet. Self-closing doors minimize human error by ensuring that the cabinet is always properly closed, even when staff are distracted or forgetful. This feature is particularly important for storing flammable liquids, as it helps contain vapors and prevent accidental releases, thereby improving overall workplace safety. The self-closing mechanism uses fusible links that melt at 165°F (74°C), causing the doors to shut automatically under fire conditions.

These cabinets are designed to meet the complex requirements of handling and storing hazardous materials, providing a reliable solution for industries such as laboratories, manufacturing, and chemical processing plants. By configuring your cabinet, you ensure a safe and compliant storage solution tailored to your workplace's unique needs.

Eagle Undercounter Safety Cabinets are backed by a ten-year warranty and are FM-approved, ensuring peace of mind and reliable performance.

Specifications

Model No	1900UC
----------	--------

Prop 65

⚠️WARNING:This product can expose you to chemicals including Lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Product Options

Color:	Red
	Yellow
	Gray
	Beige
	White
Door Style:	Self-close
	Manual Close

